

D360 Partner


Simplify secure data access and collaboration with research partners

Increasingly, biopharmaceutical research is done collaboratively by scientists across multiple organizations. Multi-organization research projects require scientists to have secure access to appropriate data, as well as easy ways to share information and ideas.

D360 Partner is an affordable add-on product for D360 for external partners that allows for seamless, secure sharing of data and analysis while advancing communication and collaboration. D360 Partner provides an alternate D360 client application for external research partners that securely connects them to relevant data from the sponsor’s main D360 instance.

- External research partners share the same data views and analysis tools as the sponsor but do not have access to unauthorized data
- Partners and sponsor can share status and ideas through D360 annotations
- Partners can be on-boarded and ramped up and down quickly and easily; partner users can be added in seconds, and data restrictions and data views can be set up as standard D360 data queries in a matter of minutes

D360 and D360 Partner Setup



External research partners can use D360 Partner to access project-relevant data views provided by the sponsor company through their existing D360 installation.

Secure Data Access with D360 Partner


Data security is of paramount importance for externalized research projects, so D360 Partner provides 3 levels of data security to prevent exposure of unauthorized data:

1. The D360 data catalog and data queries are never available within the D360 Partner client to prevent any unauthorized data access.
2. D360 Partner does not allow partners to build or adjust data queries or add additional data to datasets, preventing inadvertent data access by partners.
3. Data views are created by the sponsor and provided to external partners as simple D360 Dashboard QuickSearch widgets. The data retrieved by data view is strictly controlled by a security query, which prevents inadvertent exposure of unauthorized data by the sponsor while allowing project teams to manage data views available to partners.

Research Partner Workflow

The workflow for external partners is the same as that of the sponsor, except that building and editing of data views are omitted. Partners access data by executing D360 QuickSearch widgets provided by the sponsor that automate data formatting and visualization. The data analysis toolset and access to annotation data values are the same for the sponsor and partners, allowing cross-organizational teams to use project data in a cohesive and productive manner.

To learn more about how D360 Partner is being used in leading biopharmaceutical companies, please contact sales@certara.com.



About Certara

Certara is a leading provider of decision support technology and consulting services for optimizing drug development and improving health outcomes. Certara's solutions, which span the drug development and patient care lifecycle, help increase the probability of regulatory and commercial success by using the most scientifically advanced modeling and simulation technologies and regulatory strategies. Its clients include hundreds of global biopharmaceutical companies, leading academic institutions and key regulatory agencies.

For more information visit www.certara.com or email sales@certara.com.